



Ambalal Sarabhai Enterprises Limited

Registered Office : Shantisadan, Mirzapur Road, Ahmedabad-380001.
Telephone : +9179-25507671 / 25507073, Fax : +9179-25507483, E-mail : ase@sarabhai.co.in

Ref. No. :
Date : 07.04.2023

To,
BSE Limited,
Phiroze Jeejeebhoy Towers,
Dalal Street, Fort
Mumbai – 400 001
Scrip Code: 500009

Sub: Intimation of New Product(s) Launch by CoSara Diagnostics Private Limited, a Joint Venture of a subsidiary of Ambalal Sarabhai Enterprises Limited (“the Company”)

Dear Sir,

Pursuant to Regulation 30 read with Part B of Schedule III of the SEBI (Listing Obligations and Disclosure Requirements) Regulations, 2015, we wish to inform you that **CoSara Diagnostics Private Limited**, a Joint Venture of a subsidiary of the Company and Co-Diagnostics, Inc, Salt Lake City, Utah, USA has launched new products “**RAPID DAIGNOSTICS PRODUCT RANGE**” as per attached.

You are requested to take the same on your record.

Thanking you,


For **Ambalal Sarabhai Enterprises Limited**

DAMODAR
HARGOVINDBHAI SEJPAL




Digitally signed by DAMODAR
HARGOVINDBHAI SEJPAL
Date: 2023.04.07 16:44:56 +05'30'

Damodar H. Sejpal
Company Secretary & Compliance Officer




Encl : as above



COSARA[®]
DIAGNOSTICS
PVT. LTD



**RAPID
DIAGNOSTICS
PRODUCT RANGE**



SENSITIVE, FAST & AFFORDABLE DIAGNOSTICS


PRODUCT RANGE



SARAGENE™ HIV (1/2) Card Test

3rd Generation Kit For Differential Detection of HIV 1 & HIV 2

- Coated Ag:gp41 + gp120 for HIV 1 & gp36 for HIV 2
- Specimen Type: Serum/Plasma
- Specimen Volume: 10µl
- Result Interpretation: 10 – 15 Minutes
- Storage Temperature: 2-40°C



SARAGENE™ HBsAg (HEPATM) Card Test

Detection of HBV Surface Ag (Australian Ag)


- Specimen Type: Serum/Plasma
- Specimen Volume: 2 – 3 Drops(Dropper Provided)
- Detection Limit: 0.5ng/mL
- Result Interpretation: 20 Minutes
- Storage Temperature: 2-40°C



SARAGENE™ HCV Card Test

Detection of Ab Against HCV Antigen

- Specimen Type: Serum/Plasma
- Specimen Volume: 2 – 3 Drops (Dropper Provided)
- Result Interpretation: 15 - 20 Minutes
- Storage Temperature: 2-40°C



SARAGENE™ Syphilis Card Test

Detection of Ab against *Treponema pallidum* (Syphilis)


- Specimen Type: Serum/ Plasma
- Specimen Volume: 2 – 3 drops (Dropper Provided)
- Result Interpretation: 15 - 20 Minutes
- Storage Temperature: 2-40°C



SARAGENE™ Malaria (Pf-HRP II/Pv-pLDH) Antigen Card Test

Qualitative Detection & Dedifferentiation of Malaria Ag pLDH and HRP II


- Specimen Type: Whole Blood
- Specimen Volume: 5µl
- Detection Limit: 200 Parasites/ µl and confirm to WHO (FIND) Testing Assessment Protocol.
- Result Interpretation: 20 Minutes
- Storage Temperature: 2-40°C



SARAPLEX™ Dengue Combi Pak (NS1 Ag+IgG/IgM Ab) Card Test

Qualitative Detection of Dengue Early Ag (NS1) & Ab (IgG/IgM)

- Specimen Type: Serum / Plasma
- Specimen Volume: 2 - 3 drops (Dropper Provided)
- Result Interpretation: 15-20 Minutes
- Storage Temperature: 2-40°C



SARAGENE™ Dengue Antigen - NS1 Card Test

Qualitative Detection of Dengue Early Ag (NS1) in Human Blood

- Specimen Type: Serum/Plasma
- Specimen Volume: 2 - 3 drops (Dropper Provided)
- Result Interpretation: 15-20 Minutes
- Storage Temperature: 2-40°C



SARAGENE™ Dengue Antibody - IgG / IgM Card Test

Qualitative Detection of Dengue Ab (IgG/IgM)

- Specimen Volume: 2 - 3 drops (Dropper Provided)
- Result Interpretation: 15-20 Minutes
- Storage Temperature: 2-40°C



SARAGENE™ Typhoid IgG/IgM Card Test


Detection And Differentiation of Anti

- Salmonella Typhi (S. Typhi) IgG and IgM

- Specimen Type: Serum/ Plasma
- Specimen Volume: 2 - 3 drops (Dropper Provided)
- Result Interpretation: 15 - 20 Minutes
- Storage Temperature: 2-40°C




PRODUCT RANGE



SARAGENE™ HIV (1/2) Card Test

3rd Generation Kit For Differential Detection of HIV 1 & HIV 2


- Coated Ag:gp41 + gp120 for HIV 1 & gp36 for HIV 2
- Specimen Type: Serum/Plasma
- Specimen Volume: 10µl
- Result Interpretation: 10 – 15 Minutes
- Storage Temperature: 2-40°C



SARAGENE™ HBsAg (HEPATM) Card Test

Detection of HBV Surface Ag (Australian Ag)


- Specimen Type: Serum/Plasma
- Specimen Volume: 2 – 3 Drops(Dropper Provided)
- Detection Limit: 0.5ng/mL
- Result Interpretation: 20 Minutes
- Storage Temperature: 2-40°C



SARAGENE™ HCV Card Test

Detection of Ab Against HCV Antigen


- Specimen Type: Serum/Plasma
- Specimen Volume: 2 – 3 Drops (Dropper Provided)
- Result Interpretation: 15 - 20 Minutes
- Storage Temperature: 2-40°C




SARAGENE™ Syphilis Card Test


Detection of Ab against *Treponema pallidum* (Syphilis)

- Specimen Type: Serum/ Plasma
- Specimen Volume: 2 – 3 drops (Dropper Provided)
- Result Interpretation: 15 - 20 Minutes
- Storage Temperature: 2-40°C





Product Code	Product Name	Pack Size
	SARAGENE™	
RPHI01	HIV (1/2) Card Test	50 Test
RPHB01	HBsAg (HEPATM) Card Test	50 Test
RPHC01	HCV Card Test	50 Test
RPSY01	Syphilis Card Test	50 Test
RPMA01	Malaria (Pf-HRP II/Pv-pLDH) Antigen Card Test	50 Test
RPDE01	Dengue Antigen - NS1 Card Test	25 Test
RPTY01	Typhoid IgG/IgM Card Test	25 Test
RPDE02	Dengue Antibody - IgG / IgM Card Test	25 Test
	SARAPLEX™	
RPC001	Dengue Combi Pak (NS1 Ag+IgG/IgM Ab) Card Test	10 Test



COSARA®
DIAGNOSTICS
 PVT. LTD

Reg. Office :
 Shantisadan, Mirzapur Road,
 Ahmedabad - 380 001, India.

Factory Address : Sarabhai Campus,
 Opp. Ranoli Railway Station,
 Ranoli - 391350, Dist. Vadodara, India.

www.cosara.in

General Inquiries : info@cosara.in
 Technical Support: support@cosara.in