



SONI MEDICARE LTD.

38, Kanota Bagh, Jawahar Lal Nehru Marg, JAIPUR-302 004 INDIA
CIN No:- L51397RJ1988PLC004569
Telephone : (0141) 5163700, Fax : (0141) 2564392
E-mail : sonihospital@sonihospitals.com

Date: - 21.04.2018

To,
The Manager (Department of Corporate Affairs)
Bombay Stock Exchange Limited
Phiroze Jeejeebhoy Towers
Dalal Street, Fort
Mumbai- 400001
Email ID: corp.relations@bseindia.com

Ref: BSE Scrip Code: 539378 & SCRIP ID: SML

Sub: Compliance Certificate by RTA under Regulations 7(3) SEBI (Listing Obligations and Disclosure Requirements) Regulations, 2015

Dear Sir,

Pursuant to the requirement of Regulation 7(3) of Securities and Exchange Board of India (Listing Obligations and Disclosure Requirements) Regulations, 2015, we hereby certify that all activities in relation to both physical and electronic share transfer facility are maintained by our Share Transfer Agent viz. Sharex Dynamic (India) Pvt. Ltd, having permanent SEBI Registration No:INR000002102

Kindly take note of the above.

Thanking You

For Soni Medicare Limited
CIN: L51397RJ1988PLC004569



Aakansha Khandelwal
(Company Secretary)

COMPLIANCE CERTIFICATE

For The Half Year Ended on 31-03-2018

[Pursuant to Regulation 7(3) of the SEBI (Listing Obligations and Disclosure Requirements) Regulations, 2015]

This is to certify that pursuant to Regulation 7 (3) of SEBI (Listing Obligations and Disclosure Requirements) Regulations 2015, SONI MEDICARE LIMITED has complied with the requirements pursuant to Regulation 7 (2) of the said Regulations, by ensuring that all the activities in relation to both, Physical and Electronic, share transfer facilities i.e. Dematerialization and Rematerialization are maintained through its Registrar & Share Transfer Agent, M/s. Sharex Dynamic (India) Pvt. Ltd., which is registered with the Board (SEBI), under Registration No. INR000002102 for the Half Year ended 31-03-2018.

Yours Truly,

For SONI MEDICARE LIMITED

Aekansh



(COMPLIANCE OFFICER)

For SHAREX DYNAMIC (INDIA) PVT. LTD.

ISKumar



(AUTHORISED SIGNATORY)